


# Biosensors and Diagnostics for Pandemic Response



*PRESENTED BY*

Robert Meagher – Sandia National Laboratories

## About myself and Sandia Biosciences

Robert Meagher, Ph.D.

- Principal Member of Technical Staff
- 15 years at Sandia, starting as a postdoc
- Research interests: Microfluidics, nucleic acid amplification, virus detection

### Sandia Biosciences

- Interdisciplinary research groups located both in Livermore, CA and Albuquerque, NM
- Broad mission areas are **Biodefense and Emerging Infectious Disease** and **Bioenergy**
- Traditionally has leveraged Sandia capabilities in device & systems engineering for applications in **Biosensors, Diagnostics, and Detection**, including early foundational work in microfluidics
- The Biosciences groups collaborates broadly with other programs in Sandia, including programs in systems analysis, homeland security, global biosecurity, and other areas



## Sandia Diagnostics for COVID-19 Response



- Starting in March, Sandia as well as DOE Office of Science saw an urgent need to respond to the growing COVID-19 pandemic in the USA.
- Numerous efforts within Sandia and across the DOE Labs were rapidly initiated, leveraging multiple capabilities within the national labs (not just diagnostics).
- Two specific technologies that I have been involved with:
  - Viral RNA detection - QUASR / Smart LAMP
  - Viral protein / antibody detection - SpinDx



<https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing.html>




# “QUASR” and Smart LAMP – Rapid Viral RNA detection





- From previous research (LDRD, DTRA, and NIH funding), we had been working with isothermal diagnostics (LAMP), especially focused on arbovirus detection, and worked to advance:
  - “Direct” amplification from viruses (circumventing extraction)
  - “QUASR” detection chemistry that enabled much higher specificity as well as very bright signals that could be discerned with simple instrumentation, e.g. smartphone.
- We pivoted to SARS-CoV-2 detection, and are currently seeking commercialization opportunities.
- But we’re also interested in advancing the basic science of isothermal amplifications.

Smart Trap, Smart LAMP, QUASR for arboviruses (2013-2017)



SARS-CoV-2 Detection & Diagnostics  
(saliva, swabs, environmental)




## SpinDx – Portable instrument for multiplexed immunoassays

- Sandia began developing “SpinDx”, a portable bead-based sedimentation immunoassay, around 2009
  - Numerous applications in biodefense, emerging infectious disease, bioenergy
  - Mostly direct antigen detection, but also serology and nucleic acid detection
- SARS-CoV-2 pandemic has brought serology “test strips” into common use
  - Challenges include high variability from manufacturer to manufacturer, and unclear significance of antibody results
- Sandia SpinDx enables quantitative, multi-analyte immunoassays and even parallel immunoassay + nucleic acid assay


SpinDx prototype reader




Bead-based assay principle




SARS-CoV-2 IgG Assay



Ebola virus RNA assay





### Wearable Micro-Needle Array based sensor for bioagent exposure


- Real-time multianalyte sensing
- Proteins, RNA, small molecule biomarkers in dermal interstitial fluid
- Early warning capabilities across populations for viral infections

POC(s): Ronen Polsky ([rpolsky@sandia.gov](mailto:rpolsky@sandia.gov)),  
Phil Miller ([pruille@sandia.gov](mailto:pruille@sandia.gov))



### VOC Biomarker Discovery

- Sampling methods for breath, skin volatiles, chyme
  - State of the art laboratory equipment and expertise for biomarker discovery & ID
  - Novel modality for rapid diagnostics
- POC: Joshua Whiting ([jjwhiti@sandia.gov](mailto:jjwhiti@sandia.gov))



### Portable VOC Diagnostics technologies

- 2D Micro GC system
- Custom ion mobility spectrometer
- Custom low-SWaP electronics
- Chemical collector to improve selectivity

POC: Matthew Moorman  
([mmoorma@sandia.gov](mailto:mmoorma@sandia.gov))



**BaDx**


- No-power cartridge for culture-based anthrax detection & ID
- Self-sterilizing
- Adaptable to other pathogens

POC: Jason Harper ([jcharpe@sandia.gov](mailto:jcharpe@sandia.gov))

**Surface Acoustic Wave Biosensors**

- Sensitive, scalable, multi-mode detection (proteins, nucleic acids)
- Miniature low-power electronics
- Integrated sample prep

POC: Darren Branch ([dwbranc@sandia.gov](mailto:dwbranc@sandia.gov))



Sandia IS NOT a funding agency!

- Sandia is a federally funded research and development center (FFRDC) owned by the US Department of Energy
- (Almost) all of Sandia's work is funded by federal government (DOE and non-DOE)
- Sandia works on full cost recovery
- Sandia has limited “internal funding” - Laboratory Directed Research & Development (LDRD)
  - LDRD is a congressionally-mandated program that allows National Labs to use a small percentage of funds from federal sponsors for internally directed R&D efforts
- Sandia dedicates some portion of these LDRD funds to seed collaborative projects with Academic Alliance partners

Academic Alliance LDRD funding is not an end in itself!

- The goal is to catalyze enduring, sustainable collaborations with external sponsorship, as well as talent development



## 9 Strategizing pathways to externally funded collaborations



### University leads

- Work may be basic science, but relies on a unique Sandia capability
- Certain sponsors and mechanisms that Sandia can't be lead institution for (e.g. NSF, non-profits)

### Sandia leads

- Basic or applied science
- Good fit to national lab mission and capabilities or existing program
- Sandia may rely on university partner for unique capabilities, e.g. BSL-3, hospital collaboration

Sponsor ecosystem - some are more suitable than others for national lab involvement



# Rapid Point-of-Care Testing and Monitoring

<https://diagnostics.roche.com/us/en/products/params/electsys-anti-sars-cov-2.html>

**Jacqueline Linnes, PhD**

**Marta E. Gross Assistant Professor  
Weldon School of Biomedical Engineering  
Purdue University  
[jlinnes@purdue.edu](mailto:jlinnes@purdue.edu)**





**DIAGNOSTICS & SENSORS**  
FOR  
**FUTURE PANDEMIC RESPONSE**

**PURDUE**  
UNIVERSITY

# DIAGNOSTICS: SAMPLE-TO-RESULT DETECTION

Sample



# LABORATORY PCR-BASED TESTS




- Automated Detection
  - Multiplex samples and sample pooling
  - High throughput
- Cons:
  - Samples must be shipped to lab
  - Results not immediate





<https://helpstpauls.com/2020/03/24/covid-19/virology-laboratory-at-st-pauls-first-in-canada-to-have-new-high-throughput-covid-19-testing-capability>




# RAPID DIAGNOSTIC TESTS (ANTIGEN/ANTIBODY TESTS)



# SIGNAL ENHANCED ANTIGEN/ANTIBODY TESTS




# MOLECULAR DIAGNOSTICS




# PAPER-BASED MOLECULAR DIAGNOSTICS

Combine simplicity and portability of handheld rapid tests with sensitivity and specificity of lab-based nucleic acid diagnostics






# USER STEPS



# PAPER-BASED DIAGNOSTICS FOR RESPIRATORY INFECTIONS




Mohit Verma

COVID-19 Test: Purdue University, Raytheon BBN Technologies, a Raytheon Technologies company, PortaScience Inc., Cortex Design Inc., and LaDuca LLC



# PORTABLE SMARTPHONE PLATFORM

**Add Water**



**Start Device**



**Get Result**





OmniVis

Katherine  
Clayton




Tamara  
Kinzer-Ursem




Steven Wereley

Conflict of Interest Warning: I am a Co-Founder of OmniVis, Inc.

# VISCOSITY CHANGES WITH AMPLIFIED NUCLEIC ACID



**Low Diffusivity:**  
Pathogen Detected




**High Diffusivity:**  
No Pathogen






# SMARTPHONE DIFFUSOMETER



Camera + Computer



I: Add sample



II: Heat for 30 min



III: Add nanoparticles




IV: Deposit into chip




V: Image & analyze

# WEARABLE ELECTROPHYSIOLOGY DEVICES




Sample




# M-HEALTH SPECTROSCOPY – NONINVASIVE BLOOD ANALYSIS



Young Kim  
Natural and Biological Photonics Laboratory



# WEARABLE ELECTROPHYSIOLOGY DEVICES





# WEARABLES FOR INFECTION MONITORING



(Purdue University photo/John Underwood)



Craig Goergen

CardioVascular Imaging Research Laboratory



# Accurate, Accessible Testing for the Whole Community

Jacqueline Linnes, PhD  
[jlinnes@purdue.edu](mailto:jlinnes@purdue.edu)

17



Georgia Tech Research Institute

# Thoughts on Pathogen Detection for Emerging Threats


Michael Shannon, Ph.D.  
Mike Farrell, Ph.D.  
True Merrill, Ph.D.  
Nick Speller, Ph.D.  
Miles Paca

05NOV2020




# GTRI developed a saliva-based COVID-19 test

**GOAL:** develop a *Georgia*-based, supply chain-stable, affordable, molecular real-time reverse transcription polymerase chain reaction (RT-PCR) test intended for the qualitative detection of nucleic acid from the SARS-CoV-2 virus in saliva specimens





Omega Bio-tek  
Norcross, GA

Omni International  
Kennesaw, GA



FDA Emergency Use Authorization (pending final approval):

Georgia Tech Research Institute Saliva Collection Kit for the SARS-CoV-2 Virus  
Georgia Tech SARS-CoV-2 Virus Molecular Diagnostic Test (GT COVID-19 Saliva Test)





# Our Interests – Coupling Speed with Gold Standard Assays for Pathogen Detection




Focus across pathogen families


Attacking current gaps by creating capabilities to improve testing speed, scale and fidelity



Novel Specimen Collection



Extractionless PCR



Asymptomatic Testing at Scale



Innovation in Point of Care


# Isothermal Amplification for the Detection of SARS-CoV-2

Andrew D. Ellington

Center for Systems and Synthetic Biology

University of Texas at Austin

Sandia, November 5, 2020



Erik  
Winfree




Niles  
Pierce





Peng  
Yin

## DNA/RNA amplification LAMP




## Signal processing Toehold-mediated strand exchange computation probes








## DNA/RNA amplification LAMP




## Signal processing Toehold-mediated strand exchange computation probes






## DNA/RNA amplification LAMP




## Signal processing Toehold-mediated strand exchange computation probes




## One or multiple amplicon processing





## DNA/RNA amplification LAMP




## Signal processing Toehold-mediated strand exchange computation probes





## One or multiple amplicon processing



## Visual readout



## Detection of human fecal contamination in recreational waters



*Bacteroides* HF183 LAMP-OSD




Or-gated probe captures multiple amplicons into one signal

Use of degenerate primers captures phylogenetically wide range of Zika viruses

Bhadra et al. (2018), *Viruses* 10:714;

Bhadra et al. (2018), *PLoS Negl Trop Dis* 12:e0006671




Direct detection from mosquito guts


Cell phone image detection of individual infected mosquitoes

Results align with qPCR

- Feb 24, DIY Diagnostics Mentors and Students vote to change spring curriculum
  - New spring students implement COVID-19 CDC qPCR Dx
  - Experienced students form development team for COVID-19 LAMP-OSD Dx
- Feb 29, Existing COVID LAMP assays pulled from literature and ordered
- March 2, First successful demonstration of CDC qPCR
- March 4, screening of known primer sets begin
- March 11, OSDs designed and ordered
- March 13, campus closed, project transferred to Ellington Lab




Author	Preprint Date	Gene Target	Temp	Validation
El-Tholoth et al.	2020.02.16	ORF1ab	63	Synthesized DNA
Lamb et al. (Beaumont)	2020.02.19	ORF1a (3043-3301)	65	Synthetic Patient Samples
Yu et al.	2020.02.20	ORF1ab	63	Patient Samples
Zhang et al. (NEB)	2020.02.26	ORF1a-A	65	Patient samples
Zhang et al. (NEB)	2020.02.26	ORF1a-B	65	Patient samples
Zhang et al. (NEB)	2020.02.26	ORF1a-C	65	Patient samples
Zhang et al. (NEB)	2020.02.26	Gene N-A	65	Patient samples
Zhang et al. (NEB)	2020.02.26	Gene N-B	65	Patient samples
Yang et al.	2020.03.02	ORF1ab	63	Patient samples
Yang et al.	2020.03.02	N	63	Patient samples
Yang et al.	2020.03.02	E	63	Patient samples



**Figure 2. SARS-CoV-2 LAMP-OSD assays.** OSD fluorescence measured in real-time during LAMP amplification for NB (A), 5 primer Lamb (B), and Tholoth (C) LAMP-OSD assays are depicted as amplification curves. Presence or absence of OSD fluorescence visually observed at assay endpoint for NB (D), Lamb (E), and Tholoth (F) LAMP-OSD assays are depicted as images of reaction tubes. NB LAMP-OSD assays were seeded with indicated copies of SARS-CoV-2 N RNA or MERS-CoV N RNA or no templates. Lamb and Tholoth LAMP-OSD assays were seeded with indicated copies of gBlock DNA templates.

We also have newer primer sets that reduce detection time to 20' - 30'



**Figure 6. LAMP-OSD analysis of human saliva containing SARS-CoV-2 virions.** Indicated virion amounts were spiked in TE buffer or human saliva and added to individual or multiplex (Mx) LAMP-OSD assays. Endpoint images of OSD fluorescence are depicted for Tholoth (A), 6 primer LAMB (B), and NB (C) individual LAMP-OSD assays and Tholoth+NB Mx LAMP-OSD assays (D).

- Spike inactivated virions into saliva
- Add to LAMP-OSD assay
- That's it. Really
- Sensitivity could be better ... but keep this Figure in mind going forward




# Direct saliva analysis within 30 to 40 min


1. 3  $\mu$ L of irradiated SARS- CoV-2 virions and saliva directly added to a total of 25  $\mu$ L multiplex LAMP-OSD reactions
2. Reactions incubated at 65 °C for 30 to 40 min
3. OSD fluorescence imaged at endpoint

Mx LAMP-OSD 1 = one tube NB and Lamb6 multiplex assay  
Mx LAMP-OSD 2 = one tube NB and Tholoth multiplex assay

## Commercial RT-LAMP enzymes



## In-house engineered RT-LAMP enzymes



The addition of in-house engineered enzymes improves performance

	<b>SARS-CoV-2 DETECTR</b>	<b>SARS-CoV-2 SHERLOCK</b>	<b>CDC SARS-CoV2 qRT-PCR</b>	<b>RT-LAMP-OSD</b>
Target	N gene & E gene  (N gene gRNA compatible with CDC N2 amplicon, E gene compatible with WHO protocol)	S gene & Orf1ab gene	N-gene (3 amplicons)	Multiplex: ORF1ab, N gene
Sample control	RNase P	None	RNase P	Gapd
Limit of Detection	10-50 copies/μl input	10-100 copies/μl input	1-3.16 copies/μL input	2-20 copies/μL input (total 5 μL)
Assay reaction time	~40 min	~60 min	~120 minutes	60-90 min
Assay components	RT-LAMP (62°C, 30 min) Cas12 (37°C, 10 min) Lateral flow (RT, 2 min)	RT-RPA (42°C, 25 min) IVT + Cas13 (37°C, 30 min) Lateral flow (RT, 2 min)	UDG digestion (25°C, 2 min), reverse transcription (50°C, 15 min), denature (95°C, 2 min) amplification, (95°C, 3 sec; 55°C 30 sec; 45 cycles)	One pot RT-LAMP-OSD using either commercial enzymes Bst2.0 + RTx or using open source enzymes Bst-LF + FeRT
Heavy instrumentation required	No	No	Yes	No


Detectr, Sherlock, and qRT-PCR metrics adapted from [https://mammoth.bio/wp-content/uploads/2020/04/200423-A-protocol-for-rapid-detection-of-SARS-CoV-2-using-CRISPR-diagnostics\\_3.pdf](https://mammoth.bio/wp-content/uploads/2020/04/200423-A-protocol-for-rapid-detection-of-SARS-CoV-2-using-CRISPR-diagnostics_3.pdf)

# Cost estimate for SARS-CoV-2 test

Rough estimates based on market cost of small batch reagents


Reagent	Cost per tube (\$)	Cost (\$) per test (4 tubes)
LAMP Primers	0.11	0.44
OSD probe	0.05	0.08
Bst 2.0	0.57	2.28
RT	1.1	4.4
Betaine	0.04	0.16
dNTP	0.24	0.96
Total	2.1	8.32

- Enzymes are the most expensive ingredient.
- This cost might be reduced by using engineered enzyme alternatives



1. SARS-2 Mx RT-LAMP-OSD no primers (control)
2. SARS-2 Mx RT-LAMP-OSD with primers (test)
3. Human *gapd* no primers (control)
4. Human *gapd* with primers (test)





We have a total of 306 cases of plasmid distribution (overall)



*All thanks to Dr. Jenny Molloy at Cambridge for spearheading these efforts!*

[Home](#)[For Researchers](#)[For Partners](#)[Symposium](#)

**Reagents and  
Protocols, Share and  
Discussion**



**Forum.reclone.org**

A GLOBAL COLLABORATION FOR EQUITABLE ACCESS  
TO BIOTECHNOLOGY

# Reagent Collaboration Network

We are scientists from across the world seeking to support each other in accessing reagents for research and

**Reclone.org**

Reclone.org

## Welcome to the Reagent Collaboration Network Forum

We are scientists from across the world seeking to support each other in accessing reagents for research and diagnostics, through sharing and partnership.

Our current focus is addressing the needs of scientists researching and diagnosing COVID-19.



### Discuss with peers

Here on the Forum you can find information on protocols, available materials and discussion groups for a range of topics.



### Find and share reagents

Visit the [Shared Reagents category](#) to find what reagents are accessible within the network and how to obtain them.



### More About Us

Visit the [Reclone.org website](#) to find out about us and how to get involved. You can

 Sign up for updates

- Just about everything: Sanchita Bhadra
- Enzyme engineering: Andre Maranhao, David Walker, Inyup Paik
- Environmental screening: Dalton Towers, David Brown



Funding: NSF RAPID, NIBIB  
Supplement, Welch Foundation